

Product datasheet for **GA111574**

Human ZNF398 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF398 CRISPRa kit - CRISPR gene activation of human zinc finger protein 398
Symbol:	ZNF398
UniProt ID:	Q8TD17
Synonyms:	P51; P71; ZER6
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)