

## Product datasheet for **GA111576**

### Human TXNDC16 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                  |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                     |
| Description:  | TXNDC16 CRISPRa kit - CRISPR gene activation of human thioredoxin domain containing 16                                                                                                                                                                                                                           |
| Symbol:       | TXNDC16                                                                                                                                                                                                                                                                                                          |
| Components:   | GA111576G1, TXNDC16 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA111576G2, TXNDC16 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA111576G3, TXNDC16 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q9P2K2</a></u>                                                                                                                                                                                                                                                                                    |
| Synonyms:     | ERp90; KIAA1344                                                                                                                                                                                                                                                                                                  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                          |

#### Product images:



[View online »](#)