

Product datasheet for **GA111594**

Human ZNF319 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF319 CRISPRa kit – CRISPR gene activation of human zinc finger protein 319
Symbol:	ZNF319
Components:	GA111594G1, ZNF319 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111594G2, ZNF319 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111594G3, ZNF319 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9P2F9</u>
Synonyms:	ZFP319
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

