

## Product datasheet for **GA111608**

### Human **KCNT1** activation kit by **CRISPRa**

#### Product data:

Product Type:	CRISPRa Kits
Description:	KCNT1 CRISPRa kit - CRISPR gene activation of human potassium sodium-activated channel subfamily T member 1
Symbol:	KCNT1
UniProt ID:	<a href="#">Q5JUK3</a>
Synonyms:	bA100C15.2; EIEE14; ENFL5; KCa4.1; SLACK; Slo2.2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)