

Product datasheet for **GA111620**

Human PDZD4 activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | PDZD4 CRISPRa kit - CRISPR gene activation of human PDZ domain containing 4 |
| Symbol: | PDZD4 |
| Components: | GA111620G1, PDZD4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111620G2, PDZD4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111620G3, PDZD4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID: | <u>Q76G19</u> |
| Synonyms: | LNX5; LU1; PDZK4; PDZRN4L |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online »](#)