

Product datasheet for **GA111623**

Human WDR48 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	WDR48 CRISPRa kit - CRISPR gene activation of human WD repeat domain 48
Symbol:	WDR48
UniProt ID:	Q8TAF3
Synonyms:	P80; SPG60; UAF1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)