

Product datasheet for **GA111633**

Human KIAA1467 (FAM234B) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	FAM234B CRISPRa kit - CRISPR gene activation of human family with sequence similarity 234 member B
Symbol:	FAM234B
Components:	GA111633G1, KIAA1467 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111633G2, KIAA1467 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111633G3, KIAA1467 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	A2RU67
Synonyms:	KIAA1467
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)