

Product datasheet for **GA111636**

Human ZNF537 (TSHZ3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TSHZ3 CRISPRa kit - CRISPR gene activation of human teashirt zinc finger homeobox 3
Symbol:	TSHZ3
UniProt ID:	Q63HK5
Synonyms:	TSH3; ZNF537
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)