

Product datasheet for **GA111641**

Human LRFN1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LRFN1 CRISPRa kit - CRISPR gene activation of human leucine rich repeat and fibronectin type III domain containing 1
Symbol:	LRFN1
Components:	GA111641G1, LRFN1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111641G2, LRFN1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111641G3, LRFN1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9P244</u>
Synonyms:	SALM2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)