

## Product datasheet for **GA111642**

### Human ZFAT activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                         |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                            |
| Description:  | ZFAT CRISPRa kit - CRISPR gene activation of human zinc finger and AT-hook domain containing                                                                                                                                                                                                            |
| Symbol:       | ZFAT                                                                                                                                                                                                                                                                                                    |
| Components:   | GA111642G1, ZFAT gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA111642G2, ZFAT gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA111642G3, ZFAT gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q9P243</a></u>                                                                                                                                                                                                                                                                           |
| Synonyms:     | AITD3; ZFAT1; ZNF406                                                                                                                                                                                                                                                                                    |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                 |

#### Product images:



[View online »](#)