

Product datasheet for **GA111643**

Human KIAA1486 (NYAP2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	NYAP2 CRISPRa kit - CRISPR gene activation of human neuronal tyrosine-phosphorylated phosphoinositide-3-kinase adaptor 2
Symbol:	NYAP2
Components:	GA111643G1, KIAA1486 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GEI00074) GA111643G2, KIAA1486 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GEI00074) GA111643G3, KIAA1486 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GEI00074) 1 CRISPRa-Enhancer vector, SKU GEI00056 1 CRISPRa scramble vector, SKU GEI00077
UniProt ID:	<u>Q9P242</u>
Synonyms:	KIAA1486
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

