

Product datasheet for **GA111649**

Human p400 (EP400) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	EP400 CRISPRa kit - CRISPR gene activation of human EIA binding protein p400
Symbol:	EP400
Components:	GA111649G1, p400 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111649G2, p400 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111649G3, p400 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q96L91
Synonyms:	CAGH32; P400; TNRC12
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

