

Product datasheet for **GA111650**

Human ARHGAP23 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ARHGAP23 CRISPRa kit - CRISPR gene activation of human Rho GTPase activating protein 23
Symbol:	ARHGAP23
Components:	GA111650G1, ARHGAP23 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111650G2, ARHGAP23 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111650G3, ARHGAP23 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9P227</u>
Synonyms:	KIAA1501
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

