

Product datasheet for **GA111657**

Human **DHX37** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	DHX37 CRISPRa kit - CRISPR gene activation of human DEAH-box helicase 37
Symbol:	DHX37
UniProt ID:	<u>Q8IY37</u>
Synonyms:	DDX37
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)