

Product datasheet for **GA111663**

Human HCN3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HCN3 CRISPRa kit - CRISPR gene activation of human hyperpolarization activated cyclic nucleotide gated potassium channel 3
Symbol:	HCN3
UniProt ID:	Q9P1Z3
Synonyms:	KIAA1535; MGC131493
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)