

Product datasheet for **GA111667**

Human KIAA1543 (CAMSAP3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CAMSAP3 CRISPRa kit - CRISPR gene activation of human calmodulin regulated spectrin associated protein family member 3
Symbol:	CAMSAP3
UniProt ID:	Q9P1Y5
Synonyms:	KIAA1543; NEZHA; PPP1R80
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)