

Product datasheet for **GA111669**

Human PLEKHA4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PLEKHA4 CRISPRa kit - CRISPR gene activation of human pleckstrin homology domain containing A4
Symbol:	PLEKHA4
UniProt ID:	Q9H4M7
Synonyms:	PEPP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)