

Product datasheet for **GA111674**

Human **BEND3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	BEND3 CRISPRa kit - CRISPR gene activation of human BEN domain containing 3
Symbol:	BEND3
Components:	GA111674G1, BEND3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111674G2, BEND3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111674G3, BEND3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5T5X7</u>
Synonyms:	KIAA1553
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

