


Product datasheet for GA111686

Human TNRC6C activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TNRC6C CRISPRa kit - CRISPR gene activation of human trinucleotide repeat containing adaptor 6C
Symbol:	TNRC6C
Components:	GA111686G1, TNRC6C gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111686G2, TNRC6C gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111686G3, TNRC6C gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9HCJ0
Synonyms:	FLJ20015; KIAA1582
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)