

Product datasheet for **GA111689**

Human ZNF317 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF317 CRISPRa kit - CRISPR gene activation of human zinc finger protein 317
Symbol:	ZNF317
Components:	GA111689G1, ZNF317 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111689G2, ZNF317 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111689G3, ZNF317 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96PQ6</u>
Synonyms:	KIAA1588
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)