

## Product datasheet for **GA111705**

### Human ZNF529 activation kit by CRISPRa

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | ZNF529 CRISPRa kit - CRISPR gene activation of human zinc finger protein 529 |
| Symbol:       | ZNF529   |
| UniProt ID:   | <u><a href="#">Q6P280</a></u>  |
| Synonyms:     | KIAA1615   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)      |

#### Product images:



[View online >](#)