

Product datasheet for **GA111710**

Human PP4R4 (PPP4R4) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PPP4R4 CRISPRa kit - CRISPR gene activation of human protein phosphatase 4 regulatory subunit 4
Symbol:	PPP4R4
UniProt ID:	Q6NUP7
Synonyms:	CFAP14; KIAA1622; PP4R4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)