

Product datasheet for **GA111719**

Human Spectrin beta chain, brain 3 (SPTBN4) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SPTBN4 CRISPRa kit - CRISPR gene activation of human spectrin beta, non-erythrocytic 4
Symbol:	SPTBN4
Components:	GA111719G1, Spectrin beta chain, brain 3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111719G2, Spectrin beta chain, brain 3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111719G3, Spectrin beta chain, brain 3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9H254</u>
Synonyms:	QV; SPNB4; SPTBN3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

