

Product datasheet for GA111723

Human TRIB3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TRIB3 CRISPRa kit - CRISPR gene activation of human tribbles pseudokinase 3
Symbol:	TRIB3
UniProt ID:	Q96RU7
Synonyms:	C20orf97; NIPK; SINK; SKIP3; TRB3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)