

Product datasheet for **GA111725**

Human RB associated KRAB repressor (RBAK) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	RBAK CRISPRa kit - CRISPR gene activation of human RB associated KRAB zinc finger
Symbol:	RBAK
Components:	GA111725G1, RB associated KRAB repressor gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111725G2, RB associated KRAB repressor gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111725G3, RB associated KRAB repressor gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9NYW8
Synonyms:	ZNF769
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)