

Product datasheet for **GA111730**

Human RAB40C activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	RAB40C CRISPRa kit - CRISPR gene activation of human RAB40C, member RAS oncogene family
Symbol:	RAB40C
Components:	GA111730G1, RAB40C gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111730G2, RAB40C gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111730G3, RAB40C gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96S21</u>
Synonyms:	RARL; RASL8C
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)