

Product datasheet for **GA111760**

Human SQOR activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SQOR CRISPRa kit - CRISPR gene activation of human sulfide quinone oxidoreductase
Symbol:	SQOR
Components:	GA111760G1, SQOR gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111760G2, SQOR gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111760G3, SQOR gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9Y6N5</u>
Synonyms:	CGI-44; PRO1975; SQOR
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)