

Product datasheet for **GA111788**

Human KMT2C activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KMT2C CRISPRa kit - CRISPR gene activation of human lysine methyltransferase 2C
Symbol:	KMT2C
Components:	GA111788G1, KMT2C gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111788G2, KMT2C gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111788G3, KMT2C gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8NEZ4</u>
Synonyms:	DKFZp686C08112; FLJ12625; FLJ38309; HALR; KIAA1506
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)