

Product datasheet for **GA111792**

Human SAPAP3 (DLGAP3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DLGAP3 CRISPRa kit - CRISPR gene activation of human DLG associated protein 3
Symbol:	DLGAP3
UniProt ID:	<u>O95886</u>
Synonyms:	DAP3; SAPAP3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)