

Product datasheet for **GA111794**

Human SELK (SELENOK) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SELENOK CRISPRa kit - CRISPR gene activation of human selenoprotein K
Symbol:	SELENOK
Components:	GA111794G1, SELK gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111794G2, SELK gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111794G3, SELK gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9Y6D0</u>
Synonyms:	HSPC030; HSPC297
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)