

Product datasheet for **GA111795**

Human SINHCAF activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SINHCAF CRISPRa kit - CRISPR gene activation of human SIN3-HDAC complex associated factor
Symbol:	SINHCAF
UniProt ID:	Q9NP50
Synonyms:	C12orf14; L4; TERA
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)