

Product datasheet for **GA111797**

Human DMRT3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DMRT3 CRISPRa kit – CRISPR gene activation of human doublesex and mab-3 related transcription factor 3
Symbol:	DMRT3
Components:	GA111797G1, DMRT3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111797G2, DMRT3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111797G3, DMRT3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9NQL9</u>
Synonyms:	DMRTA3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

