

Product datasheet for **GA111808**

Human **PGAP6** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM8A CRISPRa kit - CRISPR gene activation of human transmembrane protein 8A
Symbol:	PGAP6
Components:	GA111808G1, PGAP6 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111808G2, PGAP6 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111808G3, PGAP6 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9HCN3</u>
Synonyms:	M83; TMEM6; TMEM8
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)