


Product datasheet for **GA111812**

Human HIVEP3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HIVEP3 CRISPRa kit - CRISPR gene activation of human human immunodeficiency virus type I enhancer binding protein 3
Symbol:	HIVEP3
UniProt ID:	Q5T1R4
Synonyms:	KBP-1; KBP1; KRC; Schnurri-3; SHN3; ZAS3; ZNF40C
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)