

Product datasheet for **GA111815**

Human TLNRD1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TLNRD1 CRISPRa kit – CRISPR gene activation of human talin rod domain containing 1
Symbol:	TLNRD1
Components:	GA111815G1, TLNRD1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111815G2, TLNRD1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111815G3, TLNRD1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9HIK6</u>
Synonyms:	MGC99595
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

