

Product datasheet for **GA111824**

Human PLEKHA1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PLEKHA1 CRISPRa kit - CRISPR gene activation of human pleckstrin homology domain containing A1
Symbol:	PLEKHA1
UniProt ID:	Q9HB21
Synonyms:	TAPP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)