


Product datasheet for **GA111825**

Human PLEKHA2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PLEKHA2 CRISPRa kit - CRISPR gene activation of human pleckstrin homology domain containing A2
Symbol:	PLEKHA2
Components:	GA111825G1, PLEKHA2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111825G2, PLEKHA2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111825G3, PLEKHA2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9HB19
Synonyms:	TAPP2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)