

Product datasheet for **GA111837**

Human TMEM35 (TMEM35A) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM35A CRISPRa kit - CRISPR gene activation of human transmembrane protein 35A
Symbol:	TMEM35A
Components:	GA111837G1, TMEM35 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111837G2, TMEM35 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111837G3, TMEM35 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q53FP2
Synonyms:	FLJ14084
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)