

Product datasheet for **GA111843**

Human TSKS activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TSKS CRISPRa kit - CRISPR gene activation of human testis specific serine kinase substrate
Symbol:	TSKS
Components:	GA111843G1, TSKS gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111843G2, TSKS gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111843G3, TSKS gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9UJT2</u>
Synonyms:	PPPIR161; STK22SI; TSKSI; TSSKS
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

