

Product datasheet for **GA111874**

Human **KCNK15** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	KCNK15 CRISPRa kit - CRISPR gene activation of human potassium two pore domain channel subfamily K member 15
Symbol:	KCNK15
UniProt ID:	Q9H427
Synonyms:	dJ781B1.1; K2p15.1; KCNK11; KCNK14; KT3.3; TASK-5; TASK5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)