


Product datasheet for **GA111875**

Human DHX35 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DHX35 CRISPRa kit - CRISPR gene activation of human DEAH-box helicase 35
Symbol:	DHX35
Components:	GA111875G1, DHX35 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111875G2, DHX35 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111875G3, DHX35 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9H5Z1</u>
Synonyms:	C20orf15; DDX35; KAIA0875
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)