

## Product datasheet for **GA111882**

### Human **EEFSEC** activation kit by **CRISPRa**

#### Product data:

|               |                                                                                                                                                                                                                                                                                                               |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| Description:  | EEFSEC CRISPRa kit - CRISPR gene activation of human eukaryotic elongation factor, selenocysteine-tRNA specific                                                                                                                                                                                               |
| Symbol:       | EEFSEC                                                                                                                                                                                                                                                                                                        |
| Components:   | GA111882G1, EEFSEC gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA111882G2, EEFSEC gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA111882G3, EEFSEC gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u>P57772</u>                                                                                                                                                                                                                                                                                                 |
| Synonyms:     | EFSEC; SELB                                                                                                                                                                                                                                                                                                   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:



[View online »](#)