

Product datasheet for **GA111889**

Human **SLC22A23** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	SLC22A23 CRISPRa kit - CRISPR gene activation of human solute carrier family 22 member 23
Symbol:	SLC22A23
Components:	GA111889G1, SLC22A23 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111889G2, SLC22A23 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111889G3, SLC22A23 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>A1A5C7</u>
Synonyms:	C6orf85
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)