

Product datasheet for **GA111896**

Human PKNOX2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PKNOX2 CRISPRa kit - CRISPR gene activation of human PBX/knotted 1 homeobox 2
Symbol:	PKNOX2
UniProt ID:	Q96KN3
Synonyms:	PREP2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)