

Product datasheet for **GA111902**

Human **FAM38B (PIEZO2)** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	PIEZO2 CRISPRa kit - CRISPR gene activation of human piezo type mechanosensitive ion channel component 2
Symbol:	PIEZO2
Components:	GA111902G1, FAM38B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111902G2, FAM38B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111902G3, FAM38B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9H5I5</u>
Synonyms:	C18orf30; C18orf58; DA3; DA5; FAM38B; FAM38B2; HsT748; HsT771; MWKS
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

