

Product datasheet for **GA111912**

Human Muted (BLOC1S5) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	BLOC1S5 CRISPRa kit - CRISPR gene activation of human biogenesis of lysosomal organelles complex 1 subunit 5
Symbol:	BLOC1S5
Components:	GA111912G1, Muted gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111912G2, Muted gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111912G3, Muted gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8TDH9</u>
Synonyms:	BLOS5; MU; MUTED
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)