

Product datasheet for **GA111922**

Human XPNPEP3 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | XPNPEP3 CRISPRa kit - CRISPR gene activation of human X-prolyl aminopeptidase 3 |
| Symbol: | XPNPEP3 |
| UniProt ID: | <u>Q9NQH7</u> |
| Synonyms: | APP3; ICP55; NPHPL1 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)