

Product datasheet for **GA111934**

Human DMRTB1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DMRTB1 CRISPRa kit - CRISPR gene activation of human DMRT like family B with proline rich C-terminal 1
Symbol:	DMRTB1
Components:	GA111934G1, DMRTB1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111934G2, DMRTB1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111934G3, DMRTB1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96MA1</u>
Synonyms:	DMRT-like family B with proline-rich C-terminal, 1; FLJ40434; OTTHUMP00000009751
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

