


Product datasheet for **GA111971**

Human EPB41L4A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	EPB41L4A CRISPRa kit - CRISPR gene activation of human erythrocyte membrane protein band 4.1 like 4A
Symbol:	EPB41L4A
Components:	GA111971G1, EPB41L4A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111971G2, EPB41L4A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111971G3, EPB41L4A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9HCS5</u>
Synonyms:	EPB41L4; NBL4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)