

Product datasheet for **GA111974**

Human Netrin G2 Ligand (LRRC4) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LRRC4 CRISPRa kit - CRISPR gene activation of human leucine rich repeat containing 4
Symbol:	LRRC4
Components:	GA111974G1, Netrin G2 Ligand gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111974G2, Netrin G2 Ligand gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111974G3, Netrin G2 Ligand gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9HBW1
Synonyms:	NAG14; NGL-2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

